


MT-7716 free base

Chemical Properties

CAS No. :	610323-32-5
Formula:	C27H28N4O2
Molecular Weight:	440.54
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	MT-7716 free base is a selective non-peptide agonist of nociceptin receptor (NOP). It is used for promising potential treatment drugs for alcohol abuse and relapse prevention.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2699 mL	11.3497 mL	22.6994 mL
5 mM	0.454 mL	2.2699 mL	4.5399 mL
10 mM	0.227 mL	1.135 mL	2.2699 mL
50 mM	0.0454 mL	0.227 mL	0.454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

de Guglielmo G, et al, MT-7716, a potent NOP receptor agonist, preferentially reduces ethanol seeking and reinforcement in post-dependent rats. Addict Biol. 2015 Jul;20(4):643-51.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481