

NUCC-390 dihydrochloride (1060524-97-1 free base)

Chemical Properties


CAS No. :

Formula: C₂₃H₃₅Cl₂N₅O

Molecular Weight: 468.46

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	NUCC-390 dihydrochloride is selective agonist of small-molecule CXCR4 receptor.
Targets(IC50)	CXCR
In vitro	NUCC-390 dihydrochloride (10 μ M;?pre-treatment 30 mins) leads to increased levels of pERK, it has the capability of stimulating signaling activity downstream of CXCR4 receptors[1].NUCC-390 (10 μ M;?2 hours) can induce CXCR4 receptor internalization, and non-treated cells exhibit some diffuse expression of CXCR4-YFP throughout the cytosol and clear expression in the cell membrane in HEK cells[1].
In vivo	NUCC-390 dihydrochloride?contributes to the functional and anatomical recovery of the neuromuscular junction (NMJ) following an acute nerve terminal damage by α -LTx in CD-1 mice.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1347 mL	10.6733 mL	21.3465 mL
5 mM	0.4269 mL	2.1347 mL	4.2693 mL
10 mM	0.2135 mL	1.0673 mL	2.1347 mL
50 mM	0.0427 mL	0.2135 mL	0.4269 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Mishra RK, et al. Discovery and characterization of novel small-molecule CXCR4 receptor agonists and antagonists. Sci Rep. 2016 Jul 26;6:30155.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481