


Capnolactone

Chemical Properties

CAS No. : 53771-63-4
Formula: C₁₉H₁₈O₅
Molecular Weight: 326.348
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Capnolactone is a useful organic compound for research related to life sciences. The catalog number is T124692 and the CAS number is 53771-63-4.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0642 mL	15.321 mL	30.6419 mL
5 mM	0.6128 mL	3.0642 mL	6.1284 mL
10 mM	0.3064 mL	1.5321 mL	3.0642 mL
50 mM	0.0613 mL	0.3064 mL	0.6128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481