


Guaiaretic acid

Chemical Properties

CAS No. : 500-40-3
Formula: C₂₀H₂₄O₄
Molecular Weight: 328.408
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Guaiaretic acid is a useful organic compound for research related to life sciences. The catalog number is T124764 and the CAS number is 500-40-3.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.045 mL	15.2249 mL	30.4497 mL
5 mM	0.609 mL	3.045 mL	6.0899 mL
10 mM	0.3045 mL	1.5225 mL	3.045 mL
50 mM	0.0609 mL	0.3045 mL	0.609 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481