

Piscidic acid

Chemical Properties


CAS No. :

Formula: C₁₁H₁₂O₇

Molecular Weight: 256.21

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Piscidic acid is a useful organic compound for research related to life sciences and the catalog number is T124820.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.903 mL	19.5152 mL	39.0305 mL
5 mM	0.7806 mL	3.903 mL	7.8061 mL
10 mM	0.3903 mL	1.9515 mL	3.903 mL
50 mM	0.0781 mL	0.3903 mL	0.7806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481