Data Sheet (Cat.No.T125148)



trans-Coutaric acid

Chemical Properties


CAS No.:

Formula: C13H12O8

Molecular Weight: 296.231

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	trans-Coutaric acid is a useful organic compound for research related to life sciences
	and the catalog number is T125148.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3758 mL	16.8788 mL	33.7576 mL
5 mM	0.6752 mL	3.3758 mL	6.7515 mL
10 mM	0.3376 mL	1.6879 m ∟	3.3758 mL
50 mM	0.0675 mL	0.3376 mL	0.6752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com