


Methyl arachidonate

Chemical Properties

CAS No. : 2566-89-4
Formula: C₂₁H₃₄O₂
Molecular Weight: 318.501
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Methyl arachidonate is a useful organic compound for research related to life sciences. The catalog number is T125187 and the CAS number is 2566-89-4.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1397 mL	15.6986 mL	31.3972 mL
5 mM	0.6279 mL	3.1397 mL	6.2794 mL
10 mM	0.314 mL	1.5699 mL	3.1397 mL
50 mM	0.0628 mL	0.314 mL	0.6279 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481