

Hyperwightin E

Chemical Properties


CAS No. :

Formula: C₂₉H₂₈O₁₂

Molecular Weight: 568.531

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7589 mL	8.7946 mL	17.5892 mL
5 mM	0.3518 mL	1.7589 mL	3.5178 mL
10 mM	0.1759 mL	0.8795 mL	1.7589 mL
50 mM	0.0352 mL	0.1759 mL	0.3518 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481