

Hyperwightin B

Chemical Properties


CAS No. :

Formula: C₄₀H₃₀O₁₂

Molecular Weight: 702.668

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4231 mL	7.1157 mL	14.2314 mL
5 mM	0.2846 mL	1.4231 mL	2.8463 mL
10 mM	0.1423 mL	0.7116 mL	1.4231 mL
50 mM	0.0285 mL	0.1423 mL	0.2846 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481