


1-Monolinolenin

Chemical Properties

CAS No. : 18465-99-1
Formula: C₂₁H₃₆O₄
Molecular Weight: 352.515
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8367 mL	14.1836 mL	28.3672 mL
5 mM	0.5673 mL	2.8367 mL	5.6734 mL
10 mM	0.2837 mL	1.4184 mL	2.8367 mL
50 mM	0.0567 mL	0.2837 mL	0.5673 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481