


PSI-6206 13C,d3

Chemical Properties

CAS No. : 1256490-42-2
Formula: C₁₀H₁₃FN₂O₅
Molecular Weight: 264.23
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PSI-6206 13CD ₃ is the deuterium labeled PSI-6206. PSI-6206 is a potent and selective HCV NS5B polymerase inhibitor.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7846 mL	18.9229 mL	37.8458 mL
5 mM	0.7569 mL	3.7846 mL	7.5692 mL
10 mM	0.3785 mL	1.8923 mL	3.7846 mL
50 mM	0.0757 mL	0.3785 mL	0.7569 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481