

Primverin

Chemical Properties


CAS No. :

Formula: C₂₀H₂₈O₁₃

Molecular Weight: 476.431

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0989 mL	10.4947 mL	20.9894 mL
5 mM	0.4198 mL	2.0989 mL	4.1979 mL
10 mM	0.2099 mL	1.0495 mL	2.0989 mL
50 mM	0.042 mL	0.2099 mL	0.4198 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481