


Stachybotrin F

Chemical Properties

CAS No. : 1510813-73-6
Formula: C₂₉H₃₉NO₈
Molecular Weight: 529.63
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Stachybotrin F is a useful organic compound for research related to life sciences. The catalog number is T125774 and the CAS number is 1510813-73-6.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8881 mL	9.4406 mL	18.8811 mL
5 mM	0.3776 mL	1.8881 mL	3.7762 mL
10 mM	0.1888 mL	0.9441 mL	1.8881 mL
50 mM	0.0378 mL	0.1888 mL	0.3776 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481