

Ginsenoyne K

Chemical Properties


CAS No. :

Formula: C₁₇H₂₄O₃

Molecular Weight: 276.376

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ginsenoyne K is a useful organic compound for research related to life sciences and the catalog number is T125902.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6182 mL	18.091 mL	36.1821 mL
5 mM	0.7236 mL	3.6182 mL	7.2364 mL
10 mM	0.3618 mL	1.8091 mL	3.6182 mL
50 mM	0.0724 mL	0.3618 mL	0.7236 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481