

Ajugamarin F4

Chemical Properties


CAS No. :

Formula: C₂₉H₄₂O₉

Molecular Weight: 534.646

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8704 mL	9.3519 mL	18.7038 mL
5 mM	0.3741 mL	1.8704 mL	3.7408 mL
10 mM	0.187 mL	0.9352 mL	1.8704 mL
50 mM	0.0374 mL	0.187 mL	0.3741 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481