


Syringinocide

Chemical Properties

CAS No. : 115124-95-3
Formula: C₂₃H₃₄O₁₄
Molecular Weight: 534.511
Appearance: no data available
Storage: store at low temperature
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Syringinocide is a useful organic compound for research related to life sciences. The catalog number is T126330 and the CAS number is 115124-95-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8709 mL	9.3544 mL	18.7087 mL
5 mM	0.3742 mL	1.8709 mL	3.7417 mL
10 mM	0.1871 mL	0.9354 mL	1.8709 mL
50 mM	0.0374 mL	0.1871 mL	0.3742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481