Data Sheet (Cat.No.T1271)



Pronetalol

Chemical Properties

CAS No.: 54-80-8

Formula: C15H19NO

Molecular Weight: 229.32

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Pronethalol ((±)-Pronethalol), a non-selective β-adrenergic antagonist, effectively inhibits Sox2 expression and offers protection against Digitalis-induced ventricular arrhythmias. Additionally, it can reverse such arrhythmias and limit the development of cerebral arteriovenous malformations (AVMs) [1] [2].
In vitro	Pronethalol (2, 10, 20 μ M) represses EGFP expression in a dose- and time-dependent manner in ReNcell VM cells. Pronethalol (10 μ M; 2 days) reduces the expression of Sox2 to less than 10% after 2 days of treatment [2].
In vivo	Pronethalol (0.15 mg/g; daily; for 14 days) stabilizes endothelial differentiation, lumen formation which limits cerebral arteriovenous malformation (AVMs) in Mgp -/- mice [2].

Solubility Information

Solubility	DMSO: 10 mM,
	(< 1 mg/ml refers to the product slightly soluble or insoluble)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3607 mL	21.8036 mL	43.6072 mL
5 mM	0.8721 mL	4.3607 mL	8.7214 mL
10 mM	0.4361 mL	2.1804 mL	4.3607 mL
50 mM	0.0872 mL	0.4361 mL	0.8721 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Aroesty JM, et al. Am Heart J. 1966 Apr;71(4):503-508.

Page 1 of 2 www.targetmol.com



Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 2 of 2 www.targetmol.com