


Orphenadrine hydrochloride

Chemical Properties

CAS No. : 341-69-5
 Formula: C₁₈H₂₃NO·HCl
 Molecular Weight: 305.85
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Orphenadrine hydrochloride (Mephenamin) is a muscarinic antagonist used to treat drug-induced parkinsonism and to relieve pain from muscle spasm.
Targets(IC50)	Sodium Channel,NMDAR,Norepinephrine,iGluR,Histamine Receptor

Solubility Information

Solubility	DMSO: 45 mg/mL (147.13 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2696 mL	16.3479 mL	32.6958 mL
5 mM	0.6539 mL	3.2696 mL	6.5392 mL
10 mM	0.327 mL	1.6348 mL	3.2696 mL
50 mM	0.0654 mL	0.327 mL	0.6539 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kornhuber J, et al. J Neural Transm Gen Sect. 1995;102(3):237-46.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481