


TC-MCH 7c

Chemical Properties

CAS No. :	864756-35-4
Formula:	C ₂₄ H ₂₅ FN ₂ O ₃
Molecular Weight:	408.47
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	TC-MCH 7c is an oral, selectable and blood-brain barrier penetrating MCH1R antagonist with IC ₅₀ performance of 5.6 nM against hMCH1R, and TC-MCH 7c is a phenylpyridone derivative with MCH1R K _i of 3.4 nM against human. In mice, MCH1R K _i was 3.0 nM.
Targets(IC ₅₀)	Melanin-concentrating Hormone Receptor (MCHR)
In vitro	In [Ca ²⁺] _i mobilization, TC-MCH 7c has an IC ₅₀ of 9.7 μM for MCH1R [1]. TC-MCH 7c has IC ₅₀ s of 23 nM and 9.0 μM for FLIPR and hERG, respectively[2].
In vivo	In DIO mice model, TC-MCH 7c (oral; 3-30 mg/kg; once-daily for 1.5 months) shows excellent body weight reduction in a dose-dependent manner [1]. TC-MCH 7c(oral; 3 to 30mg /kg) and 30mg /kg at 2, 15 and 24 h were 5.1, 1.8 and 0.7 μM, respectively[2].

Solubility Information

Solubility	1eq. HCl: 36.8 mg/mL (90.0 mM) DMSO: 18.0 mg/mL (44.1 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4482 mL	12.2408 mL	24.4816 mL
5 mM	0.4896 mL	2.4482 mL	4.8963 mL
10 mM	0.2448 mL	1.2241 mL	2.4482 mL
50 mM	0.049 mL	0.2448 mL	0.4896 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ito M, et al. Melanin-concentrating hormone 1-receptor antagonist suppresses body weight gain correlated with high receptor occupancy levels in diet-induced obesity mice. Eur J Pharmacol. 2009 Dec 10;624(1-3):77-83.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481