

Compound N032-0003

## Chemical Properties


CAS No. :

Formula: C<sub>32</sub>H<sub>52</sub>O<sub>3</sub>

Molecular Weight: 484.765

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## Biological Description

|             |                                                                                                                          |
|-------------|--------------------------------------------------------------------------------------------------------------------------|
| Description | Compound N032-0003 is a useful organic compound for research related to life sciences and the catalog number is T131300. |
|-------------|--------------------------------------------------------------------------------------------------------------------------|

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 2.0628 mL | 10.3142 mL | 20.6283 mL |
| 5 mM  | 0.4126 mL | 2.0628 mL  | 4.1257 mL  |
| 10 mM | 0.2063 mL | 1.0314 mL  | 2.0628 mL  |
| 50 mM | 0.0413 mL | 0.2063 mL  | 0.4126 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481