

Isounonal

Chemical Properties


CAS No. :

Formula: C₁₇H₁₂O₅

Molecular Weight: 296.278

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Isounonal is a useful organic compound for research related to life sciences and the catalog number is T131497.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3752 mL	16.8759 mL	33.7519 mL
5 mM	0.675 mL	3.3752 mL	6.7504 mL
10 mM	0.3375 mL	1.6876 mL	3.3752 mL
50 mM	0.0675 mL	0.3375 mL	0.675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481