

Perforine

Chemical Properties


CAS No. :

Formula: C₁₈H₂₅NO₅

Molecular Weight: 335.4

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9815 mL	14.9076 mL	29.8151 mL
5 mM	0.5963 mL	2.9815 mL	5.963 mL
10 mM	0.2982 mL	1.4908 mL	2.9815 mL
50 mM	0.0596 mL	0.2982 mL	0.5963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481