

(3S,5S,6R,10R,13S,17S) androstane-3,6,17-triol

Chemical Properties


CAS No. :

Formula: C₁₉H₃₂O₃

Molecular Weight: 308.462

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	androstane-3,6,17-triol is a useful organic compound for research related to life sciences and the catalog number is T131624.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2419 mL	16.2096 mL	32.4191 mL
5 mM	0.6484 mL	3.2419 mL	6.4838 mL
10 mM	0.3242 mL	1.621 mL	3.2419 mL
50 mM	0.0648 mL	0.3242 mL	0.6484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481