


Vilazodone-d8

Chemical Properties

CAS No. :	1794789-93-7
Formula:	C26H27N5O2
Molecular Weight:	449.57
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Vilazodone D8 is a deuterium-labeled vilazodone. Vilazodone is a combined inhibitor of serotonin specific reuptake (SSRI) and agonist of 5-HT1A receptor partial.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2243 mL	11.1217 mL	22.2435 mL
5 mM	0.4449 mL	2.2243 mL	4.4487 mL
10 mM	0.2224 mL	1.1122 mL	2.2243 mL
50 mM	0.0445 mL	0.2224 mL	0.4449 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481