


AP1867-3-(aminoethoxy)

Chemical Properties

CAS No. :	2127390-15-0
Formula:	C ₃₈ H ₅₀ N ₂ O ₉
Molecular Weight:	678.81
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	AP1867-3-(aminoethoxy) is a synthetic ligand for FKBP and can be used in the synthesis of PROTAC FKBP12 F36V degrader.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4732 mL	7.3658 mL	14.7317 mL
5 mM	0.2946 mL	1.4732 mL	2.9463 mL
10 mM	0.1473 mL	0.7366 mL	1.4732 mL
50 mM	0.0295 mL	0.1473 mL	0.2946 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ottis P, et al. Assessing Different E3 Ligases for Small Molecule Induced Protein Ubiquitination and Degradation. ACS Chem Biol. 2017 Oct 20;12(10):2570-2578.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481