


BMVC-8C3O

Chemical Properties

CAS No. : 1301708-12-2
Formula: C42H53I3N4O3
Molecular Weight: 1042.61
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BMVC-8C3O, a DNA G-quadruplex (G4) ligand, induces the topological conversion from non-parallel to parallel forms in human telomeric DNA G4s.
Targets(IC50)	DNA/RNA Synthesis

Solubility Information

Solubility	DMSO: 45 mg/mL (43.16 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9591 mL	4.7957 mL	9.5913 mL
5 mM	0.1918 mL	0.9591 mL	1.9183 mL
10 mM	0.0959 mL	0.4796 mL	0.9591 mL
50 mM	0.0192 mL	0.0959 mL	0.1918 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Wang ZF, et al. A novel transition pathway of ligand-induced topological conversion from hybrid forms to parallel forms of human telomeric G-quadruplexes. Nucleic Acids Res. 2016 May 5;44(8):3958-68.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481