


PROTAC FAK degrader 1

Chemical Properties

CAS No. : 2301916-69-6
Formula: C47H56F3N9O8S2
Molecular Weight: 996.13
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PROTAC FAK degrader 1 is a selective and potent degrader of focal adhesion kinase (Fak) (IC50 of 6.5 nM).
Targets(IC50)	Others,PROTACs

Solubility Information

Solubility	DMSO: 200 mg/mL (200.78 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0039 mL	5.0194 mL	10.0389 mL
5 mM	0.2008 mL	1.0039 mL	2.0078 mL
10 mM	0.1004 mL	0.5019 mL	1.0039 mL
50 mM	0.0201 mL	0.1004 mL	0.2008 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Cromm PM, et al. Addressing Kinase-Independent Functions of Fak via PROTAC-Mediated Degradation. J Am Chem Soc. 2018 Dec 12;140(49):17019-17026.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481