

PROTAC BET degrader-3

Chemical Properties


CAS No. :

Formula: C53H64N12O9S

Molecular Weight: 1045.22

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PROTAC BET Degrader-3 is a potent degrader OF BET based on PROTAC.
Targets(IC50)	Epigenetic Reader Domain

Solubility Information

Solubility	DMSO: 100 mg/mL (95.67 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9567 mL	4.7837 mL	9.5674 mL
5 mM	0.1913 mL	0.9567 mL	1.9135 mL
10 mM	0.0957 mL	0.4784 mL	0.9567 mL
50 mM	0.0191 mL	0.0957 mL	0.1913 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481