


Bragsin1

Chemical Properties

CAS No. :	369631-68-5
Formula:	C ₁₁ H ₆ F ₃ NO ₄
Molecular Weight:	273.16
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Bragsin1 is a potent and selective and noncompetitive ArfGEF BRAG2 inhibitor, inhibits Arf GTPase activation(IC ₅₀ : 3 µM),and with anti-cancer activity.
Targets(IC ₅₀)	Others

Solubility Information

Solubility	DMSO: 62.5 mg/mL (228.80 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6609 mL	18.3043 mL	36.6086 mL
5 mM	0.7322 mL	3.6609 mL	7.3217 mL
10 mM	0.3661 mL	1.8304 mL	3.6609 mL
50 mM	0.0732 mL	0.3661 mL	0.7322 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Nawrotek A, et al. PH-domain-binding inhibitors of nucleotide exchange factor BRAG2 disrupt Arf GTPase signaling. Nat Chem Biol. 2019 Feb 11. doi: 10.1038/s41589-019-0228-3.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481