


Pidotimod

Chemical Properties

CAS No. : 121808-62-6
Formula: C₉H₁₂N₂O₄S
Molecular Weight: 244.27
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pidotimod is a synthetic dipeptide with immunomodulatory properties.
Targets(IC50)	Others,Antibacterial
In vitro	Pidotimod is a synthetic dipeptide molecule with biological and immunological activity on both the adaptive and the innate immune responses. In vitro studies, both from animal and human specimens, have documented a good activity on innate and adaptive immune responses and have been confirmed by in vivo studies. These activities have been applied in clinical studies demonstrating the efficacy of pidotimod in reducing the rate of recurrent infections of the upper respiratory and urinary tracts in children. The same results are obtained in recurrent respiratory tract infections in adults.

Solubility Information

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble), DMSO: 46 mg/mL (188.3 mM), H ₂ O: 10 mg/mL (40.93 mM), < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0938 mL	20.4692 mL	40.9383 mL
5 mM	0.8188 mL	4.0938 mL	8.1877 mL
10 mM	0.4094 mL	2.0469 mL	4.0938 mL
50 mM	0.0819 mL	0.4094 mL	0.8188 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Clemente E, et al. Arzneimittelforschung, 1994, 44(12A), 1490-1494.

Riboldi P, et al. Int J Immunopathol Pharmacol, 2009, 22(2), 255-262.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 36 Washington Street, Wellesley Hills, MA 02481