


GSK-626616

Chemical Properties

CAS No. : 1025821-33-3
Formula: C₁₈H₁₀Cl₂N₄O₅
Molecular Weight: 401.27
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	GSK-626616 is a potent and orally bioavailable DYRK3 inhibitor (IC ₅₀ : 0.7 nM). It inhibits other members of the DYRK family with similar potency. GSK-626616 is a potential therapy for the treatment of anemia.
Targets(IC ₅₀)	DYRK
In vitro	GSK-626616 (serum-deprived HeLa cells) reduces EGF-induced phosphorylation of PRAS40 at Thr246. GSK-626616 enhances the binding of PRAS40 to mTORC1[2]. GSK-626616 abolishes the phosphorylation of S6K1 at Thr389 in nonstimulated HeLa cells. GSK-626616 reduces the phosphorylation of S6K1 at Thr389 in EGF- and insulin-stimulated HeLa cells, showing that mTORC1 activity is impaired[2].

Solubility Information

Solubility	DMSO: 49 mg/mL (122.11 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4921 mL	12.4604 mL	24.9209 mL
5 mM	0.4984 mL	2.4921 mL	4.9842 mL
10 mM	0.2492 mL	1.246 mL	2.4921 mL
50 mM	0.0498 mL	0.2492 mL	0.4984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Connie L. Erickson-Miller, et al. GSK626616: A DYRK3 Inhibitor as a Potential New Therapy for the Treatment of Anemia. Blood 2007 110:510;

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481