Data Sheet (Cat.No.T15862)



m-PEG3-aminooxy-Boc

Chemical Properties

CAS No.: 1835759-87-9

Formula: C12H25NO6

Molecular Weight: 279.33

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	m-PEG3-ONHBoc is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs[1].
Targets(IC50)	Others
In vitro	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

90	1mg	5mg	10mg	
1 mM	3.580 mL	17.900 mL	35.7999 mL	
5 mM	0.716 mL	3.580 mL	7.160 mL	
10 mM	0.358 mL	1.790 mL	3.580 mL	
50 mM	0.0716 mL	0.358 mL	0.716 mL	

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com