


Tetrahydropyranyldiethyleneglycol

Chemical Properties

CAS No. :	2163-11-3
Formula:	C ₉ H ₁₈ O ₄
Molecular Weight:	190.24
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Tetrahydropyranyldiethyleneglycol is a polyethylene glycol (PEG)-based PROTAC linker utilized in the synthesis of PROTACs[1].
Targets(IC50)	Others
In vitro	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2565 mL	26.2826 mL	52.5652 mL
5 mM	1.0513 mL	5.2565 mL	10.513 mL
10 mM	0.5257 mL	2.6283 mL	5.2565 mL
50 mM	0.1051 mL	0.5257 mL	1.0513 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481