

EA25

## Chemical Properties

CAS No. :

Formula: N/A

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

## Biological Description

|               |                                                                                        |
|---------------|----------------------------------------------------------------------------------------|
| Description   | EA25 has been identified in the clone library of the bacterial 16S rRNA gene sequence. |
| Targets(IC50) | Others                                                                                 |

## Solubility Information

|            |                                                                                 |
|------------|---------------------------------------------------------------------------------|
| Solubility | DMSO: 10 mM,<br>(< 1 mg/ml refers to the product slightly soluble or insoluble) |
|------------|---------------------------------------------------------------------------------|

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481