


Ald-Ph-NHS ester

Chemical Properties

CAS No. : 60444-78-2
Formula: C₁₂H₉NO₅
Molecular Weight: 247.2
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ald-Ph-NHS ester is a non-cleavable linker used in antibody-drug-conjugation (ADC).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0453 mL	20.2265 mL	40.4531 mL
5 mM	0.8091 mL	4.0453 mL	8.0906 mL
10 mM	0.4045 mL	2.0227 mL	4.0453 mL
50 mM	0.0809 mL	0.4045 mL	0.8091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481