


clAP1 Ligand-Linker Conjugates 11

Chemical Properties

CAS No. : 1239943-19-1
Formula: C35H43N3O7
Molecular Weight: 617.743
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	clAP1 Ligand-Linker Conjugates 11 is a chemical compound that combines an IAP ligand specifically designed for the E3 ubiquitin ligase with a PROTAC linker. This compound, clAP1 Ligand-Linker Conjugates 11, is useful in the development of SNIPERs[1], allowing for the targeted degradation of proteins.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6188 mL	8.094 mL	16.188 mL
5 mM	0.3238 mL	1.6188 mL	3.2376 mL
10 mM	0.1619 mL	0.8094 mL	1.6188 mL
50 mM	0.0324 mL	0.1619 mL	0.3238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Itoh Y, et al. Protein knockdown using methyl bestatin-ligand hybrid molecules: design and synthesis of inducers of ubiquitination-mediated degradation of cellular retinoic acid-binding proteins. J Am Chem Soc. 2010 Apr 28;132

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481