


clAP1 Ligand-Linker Conjugates 14

Chemical Properties

CAS No. : 1351169-40-8
Formula: C₃₉H₄₉N₃O₁₀
Molecular Weight: 719.82
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	clAP1 Ligand-Linker Conjugates 14 consist of an IAP ligand that targets the E3 ubiquitin ligase and a PROTAC linker. These conjugates, known as clAP1 Ligand-Linker Conjugates 14, are used in the development of SNIPERs[1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3892 mL	6.9462 mL	13.8924 mL
5 mM	0.2778 mL	1.3892 mL	2.7785 mL
10 mM	0.1389 mL	0.6946 mL	1.3892 mL
50 mM	0.0278 mL	0.1389 mL	0.2778 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Itoh Y, et al. Design, synthesis and biological evaluation of nuclear receptor-degradation inducers. Bioorg Med Chem. 2011 Nov 15;19(22):6768-78.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481