

HO-PEG-mal (MW 5000)

Chemical Properties

CAS No. : T18005

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	HO-PEG-mal (MW 5000) is a PEG-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others
In vitro	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481