# Data Sheet (Cat.No.T18195)



#### m-PEG36-OH

### **Chemical Properties**

CAS No.: T18195

Formula: C73H148O37

Molecular Weight: 1617.93

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

## **Biological Description**

Description	m-PEG36-OH is a PEG-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others
In vitro	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## **Preparing Stock Solutions**

	1mg	5mg	10mg	
1 mM	0.6181 mL	3.0904 mL	6.1807 mL	
5 mM	0.1236 mL	0.6181 mL	1.2361 mL	
10 mM	0.0618 mL	0.309 mL	0.6181 mL	
50 mM	0.0124 mL	0.0618 mL	0.1236 mL	

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

#### Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com