

TargetMol®
YOUR TARGET MOLECULES


Chemical Properties

Formula: C₁₇H₂₈O₄

Molecular Weight: 296.4

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Description	PEG3-C4-OBn, a polyethylene glycol (PEG)-based PROTAC linker, is utilized in crafting the PROTAC SGK3 degrader-1, a potent SGK3 degrader leveraging PROTAC technology [1][2].
Targets(IC50)	Others

	1mg	5mg	10mg
1 mM	3.3738 mL	16.8691 mL	33.7382 mL
5 mM	0.6748 mL	3.3738 mL	6.7476 mL
10 mM	0.3374 mL	1.6869 mL	3.3738 mL
50 mM	0.0675 mL	0.3374 mL	0.6748 mL

Reference

Tovell H, et al. Design and Characterization of SGK3-PROTAC1, an Isoform Specific SGK3 Kinase PROTAC Degradator. ACS Chem Biol. 2019 Sep 5.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481