


Alkyne-probe 1

Chemical Properties

CAS No. : T18915
Formula: C₁₄H₂₃N₅O₂
Molecular Weight: 293.36
Appearance: no data available
Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Alkyne-probe 1 is usually used as an Alkyne-labeled fluorescent or chemical probe.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4088 mL	17.0439 mL	34.0878 mL
5 mM	0.6818 mL	3.4088 mL	6.8176 mL
10 mM	0.3409 mL	1.7044 mL	3.4088 mL
50 mM	0.0682 mL	0.3409 mL	0.6818 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481