


Cy7 DiC18

Chemical Properties

CAS No. :	100068-60-8
Formula:	C ₆₃ H ₁₀₁ N ₂
Molecular Weight:	1013.39
Appearance:	no data available
Storage:	keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Cy7 DiC18 is a near-infrared, lipophilic fluorescent dye used for staining membranes.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9868 mL	4.9339 mL	9.8679 mL
5 mM	0.1974 mL	0.9868 mL	1.9736 mL
10 mM	0.0987 mL	0.4934 mL	0.9868 mL
50 mM	0.0197 mL	0.0987 mL	0.1974 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481