


Androstenediol 17-acetate

Chemical Properties

CAS No. : 5937-72-4
Formula: C₂₁H₃₂O₃
Molecular Weight: 332.48
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Androstenediol 17-acetate is a synthetic anabolic-androgenic steroid and an androgen ester.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0077 mL	15.0385 mL	30.077 mL
5 mM	0.6015 mL	3.0077 mL	6.0154 mL
10 mM	0.3008 mL	1.5038 mL	3.0077 mL
50 mM	0.0602 mL	0.3008 mL	0.6015 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481