


Cinoctramide

Chemical Properties

CAS No. : 28598-08-5
Formula: C₁₉H₂₇NO₄
Molecular Weight: 333.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Cinoctramide is a pharmaceutical intermediate.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 125 mg/mL (374.90 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9992 mL	14.9961 mL	29.9922 mL
5 mM	0.5998 mL	2.9992 mL	5.9984 mL
10 mM	0.2999 mL	1.4996 mL	2.9992 mL
50 mM	0.060 mL	0.2999 mL	0.5998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481