


Ethyl 3-anilinocrotonate

Chemical Properties

CAS No. :	6287-35-0
Formula:	C ₁₅ H ₁₅ NO ₂
Molecular Weight:	205.26
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Ethyl 3-anilinocrotonate is an agent of biochemical.
-------------	--

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8719 mL	24.3593 mL	48.7187 mL
5 mM	0.9744 mL	4.8719 mL	9.7437 mL
10 mM	0.4872 mL	2.4359 mL	4.8719 mL
50 mM	0.0974 mL	0.4872 mL	0.9744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481