


Ioxynic acid

Chemical Properties

CAS No. :	618-76-8
Formula:	C7H4I2O3
Molecular Weight:	389.91
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Ioxynic acid is an agent of biochemical.
-------------	--

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5647 mL	12.8235 mL	25.6469 mL
5 mM	0.5129 mL	2.5647 mL	5.1294 mL
10 mM	0.2565 mL	1.2823 mL	2.5647 mL
50 mM	0.0513 mL	0.2565 mL	0.5129 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481