


NSC-672121

Chemical Properties

CAS No. : 59147-84-1
Formula: C₁₃H₁₂O₃S
Molecular Weight: 248.3
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	NSC-672121 is an inhibitor of Cdc25A. It is also a Vitamine K3 analog.
-------------	--

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0274 mL	20.1369 mL	40.2739 mL
5 mM	0.8055 mL	4.0274 mL	8.0548 mL
10 mM	0.4027 mL	2.0137 mL	4.0274 mL
50 mM	0.0805 mL	0.4027 mL	0.8055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481