


Bamethan sulfate

Chemical Properties

CAS No. : 5716-20-1
Formula: C₂₄H₄₀N₂O₈S
Molecular Weight: 516.65
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description Bamethan sulfate is a Vasodilator.

Solubility Information

Solubility DMSO: Soluble,
(< 1 mg/ml refers to the product slightly soluble or insoluble)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9355 mL	9.6777 mL	19.3555 mL
5 mM	0.3871 mL	1.9355 mL	3.8711 mL
10 mM	0.1936 mL	0.9678 mL	1.9355 mL
50 mM	0.0387 mL	0.1936 mL	0.3871 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Li T, Cao J, Li Z, Wang X, He P. Broad screening and identification of β -agonists in feed and animal body fluid and tissues using ultra-high performance liquid chromatography-quadrupole-orbitrap high resolution mass

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481