


AHR agonist 3

Chemical Properties

CAS No. :	23749-58-8
Formula:	C ₁₈ H ₁₀ N ₂ O
Molecular Weight:	270.29
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	AHR agonist 3 is an agent with therapeutic activity.
-------------	--

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6997 mL	18.4986 mL	36.9973 mL
5 mM	0.7399 mL	3.6997 mL	7.3995 mL
10 mM	0.370 mL	1.8499 mL	3.6997 mL
50 mM	0.074 mL	0.370 mL	0.7399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481