


Tetrahydroxyquinone disodium

Chemical Properties

CAS No. :	1887-02-1
Formula:	C ₆ H ₂ Na ₂ O ₆
Molecular Weight:	216.06
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Tetrahydroxyquinone disodium is a bioactive chemical.
-------------	---

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6283 mL	23.1417 mL	46.2834 mL
5 mM	0.9257 mL	4.6283 mL	9.2567 mL
10 mM	0.4628 mL	2.3142 mL	4.6283 mL
50 mM	0.0926 mL	0.4628 mL	0.9257 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481